
GAIN™ medium

GAIN™ medium

GAIN™ medium EN

Single-step culture medium for the *in vitro* culture of human embryos and gametes

STERILE A

GAIN™ medium is sterilized by sterile filtration. Document reference: FP09 I79 GAIN R01 D.2 Update: 27.03.2019

USED ABBREVIATIONS

ICSI IntraCyttoplasmatic Sperm Injection
IUI Intra-Uterine Insemination
IVF In Vitro Fertilization
OCC Oocyte-Corona-Cumulus

INTENDED USE

GAIN™ medium is a ready to use single-step cell culture medium for use with human embryos and gametes. *For professional use only.*

COMPOSITION

GAIN™ medium is a bicarbonate-buffered balanced salt solution with 10 mg/liter gentamicin and 3.5 g/liter human serum albumin.

MATERIAL INCLUDED WITH THE KIT

Product code: GAIN010
» 1x 10ml GAIN™ medium

MATERIAL NOT INCLUDED WITH THE KIT

- » Well dishes
- » CO₂ incubator (37°C - 5% CO₂)
- » LAF bench (ISO Class 5)
- » Microscope
- » Syringe (e.g. 1 ml Plastipack)
- » Catheter (for embryo transfer)

PRODUCT SPECIFICATIONS

- » Chemical composition
- » pH @ 37°C and 5% CO₂: 7.20-7.45
- » Osmolality: 270-290 mOsm/kg
- » Sterility: SAL 10⁻³
- » Endotoxins < 0.25 EU/ml
- » Mouse Embryo Assay (blastocysts after 96h): ≥ 80%
- » Use of Ph Eur or USP grade products if applicable
- » The certificate of analysis and MSDS are available upon request

PRE-USE CHECKS

- » Do not use the product if it becomes cloudy, or shows any evidence of microbial contamination.
- » Do not use the product if the seal of the container is opened or defect when the product is delivered.

STORAGE AND CONSERVATION

- » Suitable for transport or short-term storage at elevated temperatures (up to 5 days at 37°C)
- » Store at 2-8°C, do not freeze before use.
- » Keep away from (sun)light.
- » The product can be used safely up to 7 days after opening, when sterile conditions are maintained, and the products are stored at 2-8°C.
- » Do not use after expiry date.

WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens.

Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination. GAIN™ medium contains the antibiotic Gentamicin Sulphate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.


METHOD

General notes before use

- » GAIN™ medium should generally be used in an incubator at 37°C set at 5% CO₂ at normal atmospheric pressure*
- * IMPORTANT: read *"Note on optimal pH and atmospheric pressure"*
- » Before use, always incubate GAIN™ medium overnight in a CO₂ incubator at 37°C with optimal % CO₂. Make sure that screw tops of dishes or flasks are loosened during equilibration.

Use of GAIN™ medium

GAIN™ medium can be used in all of the below procedures:

- » **Gametes:** Oocytes and semen can be handled/incubated in GAIN™ medium in preparation of, or during IVF/ICSI or IUI. For sperm processing, GAIN™ medium can be used to dilute and to wash semen for centrifugation. It can also be used in combination with a density gradient (SIL-Select Plus™, FertiPro NV, Belgium) and/or swim-up, following standard procedures. Note however that GAIN™ medium does not contain HEPES so in order to maintain a physiological pH, CO₂ incubation is necessary. If a HEPES-buffered medium is preferred, we recommend to use FertiCult™ Flushing medium (FertiPro N.V., Belgium).

- » **Fertilization and embryo culture:** GAIN™ medium can be used for fertilization, and for embryo culture from day 1 to expanded blastocyst stage (of course, shorter culture times are possible).

Example: a suitable culture system is to prepare 6 cm dishes with six 25-50 µl droplets of GAIN™ medium. To maintain pH, temperature and osmolality, make sure that the droplets are completely covered with light paraffin oil (FertiCult™ Mineral Oil, FertiPro NV, Belgium). Ideally, these dishes are prepared the day before use and incubated overnight in a CO₂ incubator at 37°C and 5% CO₂. Per droplet 1-4 trimmed OCCs can be cultured and inseminated with approximately 0.1 million sperms.

GAIN™ medium is a single-step medium designed for continuous embryo culture from day 1 to expanded blastocyst stage. However, if culture medium is refreshed, it is recommended to do this on day 2 or early day 3. In this case, make sure that the fresh medium is pre-incubated before transfer of the embryos.

- » GAIN™ medium can also be used for **embryo transfer**.

Note on optimal pH and atmospheric pressure

GAIN™ medium must generally be used in an incubator at 37°C and 5% CO₂. At normal atmospheric pressure at sea level (101.3 kPa, range: 99 - 103 kPa) this will result in a medium with a pH of approximately 7.28. At lower atmospheric pressure or higher altitudes, the percentage CO₂ in the incubator must be increased. Below is a table with estimated CO₂ setpoint at different air pressures. The only way to obtain certainty about the pH is to measure directly under culture conditions with a suitable well-calibrated pH meter. For optimal embryo culture the pH should be: 7.28 ± 0.05. GAIN™ medium contains a low concentration of Phenol Red (approx. 0.001 mM) which is sufficient to indicate the pH. At the correct pH the colour is a pale bright red. If the medium turns orange or yellowish the pH is too low, if it turns pink or purplish the pH is too high.

Altitude (meters)	Atmospheric air pressure	CO ₂ setpoint*
0 (sea level)	101.3 kPa	5.0%
0 - 400 m	100 - 96 kPa	5.0 - 5.2%
400 - 800 m	96 - 92 kPa	5.3 - 5.5%
800 - 1200 m	92 - 89 kPa	5.6 - 5.8%

* Calculated as follows: (pressure sea level/pressure altitude) x 5%.

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- Huisman GJ, Alberda AT, Leertveld RA, Verhoeff A, Zeilmaker GH, 1994. A comparison of in vitro fertilization results after embryo transfer after 2, 3, and 4 days of embryo culture. *Fertil Steril* 61, 970-72
- Scholtes MCW, Zeilmaker GH, 1996. A prospective, randomized study of embryo transfer results after 3 or 5 days of embryo culture in in vitro fertilization. *Fertil Steril* 65, 1245-48
- Rijnders PM, Jansen CAM, 1998. The predictive value of day 3 embryo-morphology for blastocyst formation and implantation rate at day 5 in IVF. *Human Reprod* 13, 2869-73
- Rijnders PM, Jansen CAM, 1999. Influence of group culture and culture volume on the formation of human blastocysts: a prospective randomised study. *Human Reprod* 14, 2333-7
- Macklon NS, Pieters MHEC, Hassan MA, Jeucken PHM, Eijkemans MJC, Fauser BQJM, 2002. A prospective randomized comparison of sequential versus monoculture systems for in-vitro human blastocyst development. *Human Reprod* 17, 2700-05
- Curfs MHJM, Cleine JH, Hondelink MN, Van Kamp AA, Kruse ME, Leertveld RA, Comparison of two types of Embryo transfer catheter. Poster presented at the 3rd International Alpha Congress (International Society of Clinical Embryologists), "ART, Science and Fiction", 9-11 September, New York 2001

TECHNICAL SUPPORT
ASSISTANCE TECHNIQUE
TECHNISCHE UNTERSTÜTZUNG
SUPPORTO TECNICO
AYUDA TÉCNICA
SUPPORTO TÉCNICO

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