
EN FR DE

# GAIN™ medium

## GAIN™ medium

Single-step culture medium for the *in vitro* culture of human embryos and gametes

STERILE A

GAIN™ medium is sterilized by sterile filtration. Document reference: FP09 I79 GAIN R01 D.2 Update: 27.03.2019

### USED ABBREVIATIONS

ICSI IntraCytoplasmatic Sperm Injection  
IUI Intra-Uterine Insemination  
IVF In Vitro Fertilization  
OCC Oocyte-Corona-Cumulus

### INTENDED USE

GAIN™ medium is a ready to use single-step cell culture medium for use with human embryos and gametes. *For professional use only.*

### COMPOSITION

GAIN™ medium is a bicarbonate-buffered balanced salt solution with 10 mg/liter gentamicin and 3.5 g/liter human serum albumin.

### MATERIAL INCLUDED WITH THE KIT

Product code: GAIN010  
» 1x 10ml GAIN™ medium

### MATERIAL NOT INCLUDED WITH THE KIT

» Well dishes  
» CO<sub>2</sub> incubator (37°C - 5% CO<sub>2</sub>)  
» LAF bench (ISO Class 5)  
» Microscope  
» Syringe (e.g. 1 ml Plastipack)  
» Catheter (for embryo transfer)

### PRODUCT SPECIFICATIONS

» Chemical composition  
» pH @ 37°C and 5% CO<sub>2</sub>: 7.20-7.45  
» Osmolality: 270-290 mOsm/kg  
» Sterility: SAL 10<sup>-3</sup>  
» Endotoxins < 0.25 EU/ml  
» Mouse Embryo Assay (blastocysts after 96h): ≥ 80%  
» Use of Ph Eur or USP grade products if applicable  
» The certificate of analysis and MSDS are available upon request

### PRE-USE CHECKS

» Do not use the product if it becomes cloudy, or shows any evidence of microbial contamination.  
» Do not use the product if the seal of the container is opened or defect when the product is delivered.

### STORAGE AND CONSERVATION

» Suitable for transport or short-term storage at elevated temperatures (up to 5 days at 37°C)  
» Store at 2-8°C, do not freeze before use.  
» Keep away from sunlight.  
» The product can be used safely up to 7 days after opening, when sterile conditions are maintained, and the products are stored at 2-8°C.  
» Do not use after expiry date.

### WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination. GAIN™ medium contains the antibiotic Gentamicin Sulphate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.


### METHOD

#### General notes before use

- » GAIN™ medium should generally be used in an incubator at 37°C set at 5% CO<sub>2</sub> at normal atmospheric pressure\*
- \* IMPORTANT: read "Note on optimal pH and atmospheric pressure"
- » Before use, always incubate GAIN™ medium overnight in a CO<sub>2</sub> incubator at 37°C with optimal % CO<sub>2</sub>. Make sure that screw tops of dishes or flasks are loosened during equilibration.

#### Use of GAIN™ medium

GAIN™ medium can be used in all of the below procedures:

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