EN FR DE IT

10% PVP in FertiCult™ Flushing medium

10% PVP in FertiCult™ Flushing medium

10% PVP in FertiCult™ Flushing medium

Viscous medium for reducing sperm motility

STERILE A

Doc. reference: FP09 I17 R01 D.5
Update: 18.08.2020

USED ABBREVIATIONS

ICSI Intra-cytoplasmic Sperm Injection
PVP Polyvinylpyrrolidone

GENERAL INFORMATION AND INTENDED USE

Intra-cytoplasmic sperm injection (ICSI) requires the capture of individual sperm cells in a glass pipette for injection into the oocyte.

This procedure is facilitated by first immobilizing the sperm. The most widely used method of immobilizing the sperm cells is by placing them in a viscous medium prior to nicking the tail to immobilize sperm completely.

10% PVP in FertiCult™ Flushing medium is a ready-to-use viscous medium that can be used in ICSI procedures. For professional use in human assisted reproduction techniques.

COMPOSITION

10% (w/v) dialysed Ph Eur grade polyvinylpyrrolidone (PVP) in FertiCult™ Flushing medium. FertiCult™ Flushing medium is an aqueous solution containing physiologic salts, HEPES, lactate, pyruvate, glucose and human albumin solution.

MATERIAL INCLUDED WITH THE KIT

PVP1
1 x 1ml 10% PVP in FertiCult™ Flushing medium

PVP0.2S
5 x 0.2ml 10% PVP in FertiCult™ Flushing medium

MATERIAL NEEDED BUT NOT PROVIDED

- » Incubator at 37°C (No CO₂)
- » Petri dishes
- » Glass pipettes
- » FertiCult™ Mineral Oil
- » FertiCult™ Flushing medium
- » LAF bench
- » Microscope

PRODUCT SPECIFICATIONS

- » Chemical composition
- » pH: 7,20-7,90 (Release criterion: 7,20-7,60)
- » Osmolality: 280-320 mOsm/kg
- » Viscosity: > 170 cP at 25°C
- » Sterility: sterile (SAL 10⁻³)
- » Endotoxins: < 1 EU/ml
- » Mouse Embryo Assay (blastocysts after 96h): ≥ 80% after 10 min exposure
- » Use of Ph Eur or USP grade products if applicable
- » A certificate of analysis and MSDS are available upon request

10% PVP in FertiCult™ Flushing medium has been evaluated by an independent test laboratory using the ISO 10993-5/ USP30-NF25 Agar Diffusion test for cytotoxicity. The product was considered non-cytotoxic.

Symbols as defined in ISO 15223-1:2016
Symboles tels que définis dans ISO 15223-1:2016
Symbole nach der Definition von ISO 15223-1:2016
Simboli di cui alla norma ISO 15223-1:2016
Símbolos definidos en la norma ISO 15223-1:2016
Símbolos conforme definido na norma ISO 15223-1:2016
Symbolen volgens ISO 15223-1:2016
Σύμβολα όπως ορίζονται στο ISO 15223-1:2016

PRE-USE CHECKS

- » Do not use the product if it becomes cloudy or shows any evidence of microbial contamination.
- » Do not use the product if seal of the container is opened or defect when the product is delivered.

STORAGE INSTRUCTIONS

- » Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).
- » Store between 2-8°C
- » Do not freeze before use
- » Keep away from sunlight
- » The product can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C
- » Do not use after expiry date

WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection / removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (ISO 5 environment, e.g LAF-bench) to avoid possible contamination.

10% PVP in FertiCult™ Flushing medium does not contain any antibiotics. The long term safety of ICSI, including the use of PVP, on children born following this procedure is unknown.

METHOD

Evaporation of the PVP solution on the slide or dish and the presence of human albumin may occasionally cause a small amount of precipitation in the medium.

Standard procedure

1. Warm the PVP solution to 37°C.
2. Place a small drop of PVP solution (5µl – 10µl) in a dish and cover with FertiCult™ Mineral oil.
3. Add a small volume (1µl – 2µl) of washed sperm into the centre of the PVP droplet.
4. Incubate for a (few) minute(s) to allow sperm to migrate to the periphery of the droplet.
5. Select and recover the spermatozoa for injection.

Alternative procedure with extra washing step

1. Warm the PVP solution and FertiCult™ Flushing medium to 37°C.
2. Place a small drop of PVP solution (5µl – 10µl) and 1 or more small drops FertiCult™ Flushing medium in a dish and cover with FertiCult™ Mineral Oil.
3. Add a small volume (1µl – 2µl) of washed sperm into the centre of the PVP droplet.
4. Incubate for a (few) minute(s) to allow sperm to migrate to the periphery of the droplet.
5. Select the spermatozoa for injection and nick (break) the tail of the spermatozoon with the tip of the glass pipette.
6. Transfer the spermatozoon into one of the FertiCult™ Flushing medium droplets and wash by transferring the sperm cell in and out of the flushing medium several times.
7. Suck the sperm cell into the glass pipette and use for ICSI procedure.

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- Gianaroli L., et al, ESHRE guidelines for good practice in IVF laboratories. Human Reproduction, 2000; 15 (No. 10); pp2241-2246.
- Textbook of in vitro fertilization and assisted reproduction, The Bourn Hall Guide to Clinical and Laboratory Practice 3rd Edition, Ed: Brnnsden PR, Taylor and Francis, 2005, pp. 342-343.
- Noorashikin M et al. Affordable ART for developing countries: a cost benefit comparison of low dose stimulation versus high dose GnRH antagonist protocol. Journal of Assisted Reproduction and Genetics, 2008; 25(7):297-303.
- Zainul Rashid MR, et al. GnRH Agonist and GnRH antagonist in intracytoplasmic injection cycles. Med J Malaysia, 2008, 63(2).

TECHNICAL SUPPORT ASSISTANCE TECHNIQUE TECHNISCHE UNTERSTÜTZUNG SUPPORTO TECNICO ASISTENCIA TÉCNICA SUPORTE TÉCNICO TECHNISCHE ONDERSTEUNING ΤΕΧΝΙΚΗ ΥΠΟΣΤΗΡΙΞΗ

FertiPro N.V.
Industriepark Noord 32
8730 Beernem – Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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