

IUI PLUS

UTERINE INSEMINATOR FOR WASHED SPERMATOZOA WITH MANDREL

The IUI Plus is a set for difficult cervical passage, consisting of two parts:

- A) Loading cannula
- B) Mandrel


FEATURES

- positioner for depth adjustment
- a combination of a soft, flexible intrauterine cannula and a solid cervix cannula
- markings on outer tube for easy positioning
- mandrel to prevent mucus to enter the loading cannula
- handle with luer lock for a more secure fitting to syringes

ORIGINAL PICTURE

GYNO INFO IUI PLUS

Length	A) 216 mm (± 3 mm) B) 240 mm (± 2 mm)
Tube OD	A) 1.60 mm (± 0.10 mm) B) 1.00 mm (± 0.05 mm)
Total volume	A) 0.42 ml
Total volume minus hub	A) 0.258 ml
Sterile	sterilized by gamma irradiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4220 Plus ②



GMP Belgium IUI Plus 06/2014



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS®
Practise better Medicine