

IUI

UTERINE INSEMINATOR FOR WASHED SPERMATOZOA

The IUI is an intrauterine insemination cannula with two lateral distal ports for optimal sperm distribution.

FEATURES

- positioner to prevent damage to fundus
- two lateral ports at distal end for optimal sperm distribution
- arrows on handle to indicate lateral ports at distal end
- markings on the outer tube
- handle with luer lock for a more secure fitting to syringes



ORIGINAL PICTURE

GYNO INFO IUI

Length	206 mm ($\pm$ 2 mm)
Tube OD	1.85 mm ($\pm$ 0.10 mm)
Sterile	sterilized by gamma irradiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	#4220 



GMP Belgium IUI 09/2014

