


# SUPERIOR FAS DOUBLE LUMEN FLUSH

## CUSTOM-MADE FOLLICLE ASPIRATION SET

The Superior FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.

### FEATURES

- flushing connector enables flushing and aspirating during complicated oocyte retrieval or when few mature oocytes occur
- double lumen allows simultaneous flushing and aspirating
- dual laser and embossed echogenic marked tip available which offer stable and excellent vision during ultrasound
- translucent tubing in the same size as needle for safe oocyte retrieval
- various needle sizes for ovaries with complicated and/or limited access
- arrow marking on handle for bevel point indication
- handle shape facilitates rotation
- custom-made



ORIGINAL PICTURE

### GYNO INFO SUPERIOR FAS DOUBLE LUMEN FLUSH

Sterile	sterilized by gamma irradiation
Packaging	50/shipping box (standard configuration)
Quality control	MEA and LAL
Order number	<b>#4551</b> ② (complete the configurator and receive your order number)

