

SEMTRAC C

SOFT EMBRYO TRANSFER CATHETER

The Semtrac C is a soft embryo transfer catheter.

FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- one piece economical catheter
- markings to indicate depth
- metal tube inside for better control
- handle with luer lock for a more secure fitting to syringes

ORIGINAL PICTURE

GYNO INFO SEMTRAC C

Length	252 mm ($\pm$ 2 mm)
Tube OD	1.70 mm ($\pm$ 0.05 mm)
Total volume	0.20 ml
Total volume minus hub	0.065 ml
Sterile	ETO sterilized
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#2000 C ②



GMP Belgium Semtrac C 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNÉTICS®
Practise better Medicine