

EMTRAC SET

EMBRYO TRANSFER CATHETER WITH POSITIONER

The Emtrac Set is a set for difficult cervical passage, consisting of three parts:

- A) Guiding catheter
- B) Mandrel
- C) Loading catheter


FEATURES

- positioner for depth adjustment
- a combination of a soft, flexible intrauterine catheter and a solid cervix catheter
- markings on outer tube for easy positioning
- soft loading catheter to minimize trauma
- handle with luer lock for a more secure fitting to syringes

ORIGINAL PICTURE

GYNO INFO EMTRAC SET

Length	A) 216 mm (± 3 mm) B) 240 mm (± 2 mm) C) 266 mm (± 2 mm)
Tube OD	A) 1.60 mm (± 0.10 mm) B) 1.00 mm (± 0.05 mm) C) 1.00 mm (± 0.05 mm)
Total volume	C) 0.20 ml
Total volume minus hub	C) 0.069 ml
Sterile	sterilized by radiation
Packaging	75/shipping box
Quality control	MEA and LAL
Order number	#4219 Set ②



GMP Belgium Emtrac Set 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine