

SEMTRAC SET 4.2

SOFT EMBRYO TRANSFER CATHETER

The Semtrac Set is a soft embryo transfer catheter set for difficult embryo transfers.

FEATURES

- soft bullet shaped tip to minimize trauma
- individual microscopic quality inspection
- three piece catheter for difficult transfers
- 5 **free** mandrels, to preshape the guiding catheter, in every box of 25 pieces
- reinforced shaft of loading catheter for optimal control
- handle with luer lock for a more secure fitting to syringes

ORIGINAL PICTURE

GYNO INFO SEMTRAC SET 4.2

Guiding catheter (A) length 158 mm (± 2 mm)
tube OD 2.55 mm (± 0.05 mm)


Loading catheter (C) length 232 mm (± 2 mm)
tube OD 1.70 mm (± 0.05 mm)
total volume 0.23 ml
total volume minus hub 0.081 ml

Sterile ETO sterilized

Packaging 75/shipping box

Quality control MEA and LAL

Order number **#2000 Set 4.2** ②



GMP Belgium Semtrac Set 4.2 04/2011



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine