

# Oosafe® IVF TESTED PLASTIC WARE FOR IVD USE



**CE** Oosafe® Plasticware is manufactured  
in an EN ISO 13485:2016 certified facility.

# Oosafe® IVF TESTED PLASTIC WARE FOR IVD USE

## TOP QUALITY PLASTICWARE USED BY SCIENTISTS ALL OVER THE WORLD

- Designed for wide range of applications
- Made of Medical Grade USP VI Raw Materials:
  - Polystyrene (PS)
  - Polypropylene (PP)
  - Polyethylene (PE)
- Surface treated options for a consistent hydrophilic surface
- Packaged in peel or tear-open packaging for user convenience
- Excellent thermal transfer on heated stage
- Perfect fit on benchtop incubators with short rim on the bottom (except ICSI/IMSI Dish which has totally flat bottom)
- Easy Grip and Label Area Design
- Sterilized with Gamma Irradiation
- CE marked according to 98/79/EC Directive on *In Vitro* Diagnostic Medical Devices
- Tested for each LOT:
  - One cell stage Mouse Embryo Assay (MEA > 80% Blastocyst)
  - Non-pyrogenic, Endotoxin (LAL < 0.03125 EU/device)
  - Human Sperm Survival and Motility (SMI ≥ 0.75%)
- Vacuum packed to prevent scratches



**Oosafe® 4 Well Dish**  
NON-TREATED SURFACE

MULTIPACK PS

4 pcs/pack  
120 pcs/case



Product code:  
**OOPW-FW04**

Cut corner design ensures proper orientation  
and prevents mix-ups while handling  
Well volume 1.8 ml  
Volume between wells 7 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



**Oosafe® 4 Well Dish**  
TREATED SURFACE

MULTIPACK PS

4 pcs/pack  
120 pcs/case



Product code:  
**OOPW-FW03**

Cut corner design ensures proper orientation  
and prevents mix-ups while handling  
Well volume 1.8 ml  
Volume between wells 7 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



**Oosafe® 6 Well Dish**  
with Straw Holder

SINGLEPACK PS

1 pc/pack  
240 pcs/case



Product code:  
**OOPW-SW02**

Conical flat bottom. Hands-free design: the tip of the straw  
and the bottom of the well are at the same focus level  
Well volume 0.3 ml  
Volume between wells 3 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



**Oosafe® 6 Well Dish**  
with Straw Holder

MULTIPACK PS

4 pcs/pack  
120 pcs/case



Product code:  
**OOPW-SW03**

Conical flat bottom. Hands-free design: the tip of the straw  
and the bottom of the well are at the same focus level  
Well volume 0.3 ml  
Volume between wells 3 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® Center Well Dish with Two Compartments

LABEL AREA GRIP  
SINGLEPACK PS

1 pc/pack  
240 pcs/case



Product code:  
**OOPW-CW04**

Multi-purpose dish for cell handling applications  
Center well volume 3 ml  
Each compartment volume 5 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® Center Well Dish with Two Compartments

LABEL AREA GRIP  
MULTIPACK PS

10 pcs/pack  
500 pcs/case



Product code:  
**OOPW-CW05**

Multi-purpose dish for cell handling applications  
Center well volume 3 ml  
Each compartment volume 5 ml

### Oosafe® ICSI/IMSI Dish for Sperm Selection

LABEL AREA GRIP  
MULTIPACK PS

10 pcs/pack  
500 pcs/case



Product code:  
**OOPW-IC03**

Complete flat and extra thin bottom  
for perfect image quality and thermal transfer  
Volume 12 ml

### Oosafe® ICSI/IMSI Dish for Sperm Selection

LABEL AREA GRIP  
SINGLEPACK PS

1 pc/pack  
240 pcs/case



Product code:  
**OOPW-IC02**

Complete flat and extra thin bottom  
for perfect image quality and thermal transfer  
Volume 12 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® Center Well Dish

LABEL AREA GRIP  
MULTIPACK PS

10 pcs/pack  
500 pcs/case



Product code:  
**OOPW-CW08**

New design with one outer chamber  
Center well volume 3 ml  
Outer compartment volume 10 ml

### Oosafe® 60 mm Dish,

LABEL AREA GRIP  
MULTIPACK PS

10 pcs/pack  
500 pcs/case



Product code:  
**OOPW-ST03**

Multi-purpose dish for cell culture  
and microscopy applications  
Volume 25 ml

### Oosafe® 35 mm Dish, High Wall

LABEL AREA GRIP  
MULTIPACK PS

10 pcs/pack  
500 pcs/case



Product code:  
**OOPW-TF03**

Multi-purpose dish for cell culture and microscopy  
applications, high wall with enough space  
for labels between the grips  
Volume 8 ml

### Oosafe® Synergy Dish

TREATED SURFACE  
MULTIPACK PS

4 pcs/pack  
120 pcs/case



Product code:  
**OOPW-SD04**

Unique design with micro-wells for group culture  
Volume of wash well 56 µl  
Volume of group culture well 160 µl  
Total volume 15 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® 100 mm Dish

MULTIPACK PS

10 pcs/pack  
250 pcs/case



Product code:  
**OOPW-HD10**

Multi-purpose dish for cell culture and microscopy applications, perfect for higher volume needs  
Volume 105 ml

### Oosafe® Andrology Tube, 5 mL

MULTIPACK PS

10 pcs/pack  
500 pcs/case



Product code:  
**OOPW-AT10**

Push-on cap prevents sample loss with ease off handling

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® Centrifuge Tube, 15 mL

SINGLEPACK PP

1 pc/pack  
300 pcs/case



Product code:  
**OOPW-CT01**

Suited for applications requiring greater thermal and chemical stability

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® OPU Tube, 14 mL

MULTIPACK PS

10 pcs/pack  
500 pcs/case



Product code:  
**OOPW-OT10**

Push-on cap prevents sample loss with ease off handling

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® Sperm Collection Cup, 80 mL

SINGLEPACK PP

1 pc/pack  
70 pcs/case



Product code:  
**OOPW-SC01**

Wide opening for easy sampling  
Volume 80 ml

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® PP Centrifuge Tube, 50 mL

MULTIPACK PP

25 pcs/pack  
500 pcs/case



Product code:  
**OOPW-CT25**

Ideal for higher volume centrifugation

### Oosafe® PS Centrifuge Tube, 15 mL

MULTIPACK PS

5 pcs/pack  
500 pcs/case



Product code:  
**OOPW-CT05**

Suited for procedures requiring high optical clarity

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® Flask, 50 mL

TREATED SURFACE  
MULTIPACK PS

10 pcs/pack  
160 pcs/case



Product code:  
**OOPW-FL10**

Bent neck allows to access all corners of the flask

### Oosafe® Pasteur Pipette, 1 mL

SINGLEPACK PE

1 pc/pack  
500 pcs/case



Product code:  
**OOPW-PP01**

*For quick and safe transfer of media and samples*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® Pasteur Pipette, 3 mL

SINGLEPACK PE

1 pc/pack  
500 pcs/case



Product code:  
**OOPW-PP03**

*For quick and safe transfer of media and samples*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® 1 mL in 1/100 mL Serological Pipette

SINGLEPACK PS

1 pc/pack  
500 pcs/case



Product code:  
**OOPW-SP01**

*Negative graduations for extra capacity*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® 2 mL in 1/100 mL Serological Pipette

SINGLEPACK PS

1 pc/pack  
500 pcs/case



Product code:  
**OOPW-SP02**

*Negative graduations for extra capacity*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® 5 mL in 1/10 mL Serological Pipette

SINGLEPACK PS

1 pc/pack  
250 pcs/case



Product code:  
**OOPW-SP05**

*Negative graduations for extra capacity  
and reverse graduations for user convenience*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

### Oosafe® 10 mL in 1/10 mL Serological Pipette

SINGLEPACK PS

1 pc/pack  
200 pcs/case



Product code:  
**OOPW-SP10**

*Negative graduations for extra capacity  
and reverse graduations for user convenience*

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

#### PRODUCTS AVAILABLE UPON REQUEST \*



**Oosafe® Center Well Dish  
with Two Compartments**  
ALL AROUND GRIP

MULTIPACK PS

Product code:  
**OOPW-CW03**  
10 pcs/pack  
500 pcs/case



**Oosafe® 60 mm Dish,**  
LABEL AREA GRIP  
TREATED SURFACE

MULTIPACK PS

Product code:  
**OOPW-ST05**  
10 pcs/pack  
500 pcs/case



**Oosafe® 35 mm Dish,**  
High Wall  
LABEL AREA GRIP  
TREATED SURFACE

MULTIPACK PS

Product code:  
**OOPW-TF05**  
10 pcs/pack  
500 pcs/case

\*Minimum purchase quantity of 1000 cases.

# Oosafe® IVF TESTED PLASTIC WARE FOR IVD USE



## ORDER INFORMATION

Product code	Product name	Packaging
OOPW-FW04	Oosafe® 4 Well Dish, Non-Treated Surface	MULTIPACK
OOPW-FW03	Oosafe® 4 Well Dish, Treated Surface	MULTIPACK
OOPW-SW02	Oosafe® 6 Well Dish with Straw Holder	SINGLEPACK
OOPW-SW03	Oosafe® 6 Well Dish with Straw Holder	MULTIPACK
OOPW-CW04	Oosafe® Center Well Dish with two Compartments, Label Area Grip	SINGLEPACK
OOPW-CW05	Oosafe® Center Well Dish with two Compartments, Label Area Grip	MULTIPACK
OOPW-CW08	Oosafe® Center Well Dish	MULTIPACK
OOPW-IC02	Oosafe® ICSI/IMSI Dish for Sperm Selection, Label Area Grip	SINGLEPACK
OOPW-IC03	Oosafe® ICSI/IMSI Dish for Sperm Selection, Label Area Grip	MULTIPACK
OOPW-ST03	Oosafe® 60mm Dish, Label Area Grip	MULTIPACK
OOPW-TF03	Oosafe® 35mm Dish, High Wall	MULTIPACK
OOPW-HD10	Oosafe® 100mm Dish	MULTIPACK
OOPW-SD04	Oosafe® Synergy Dish, Treated	MULTIPACK
OOPW-SC01	Oosafe® Sperm Collection Cup, 80 mL	SINGLEPACK
OOPW-CT05	Oosafe® PS Centrifuge Tube, 15 mL	MULTIPACK
OOPW-CT01	Oosafe® Centrifuge Tube, 15 mL	SINGLEPACK
OOPW-OT10	Oosafe® OPU Tube, 14 mL	MULTIPACK
OOPW-AT10	Oosafe® Andrology Tube, 5mL	MULTIPACK
OOPW-CT25	Oosafe® PP Centrifuge Tube, 50 mL	MULTIPACK
OOPW-FL10	Oosafe® Flask, 50 mL, Treated	MULTIPACK
OOPW-PP01	Oosafe® Pasteur Pipette, 1 mL	SINGLEPACK
OOPW-PP03	Oosafe® Pasteur Pipette, 3 mL	SINGLEPACK
OOPW-SP01	Oosafe® 1 mL in 1/100 mL Serological Pipette	SINGLEPACK
OOPW-SP02	Oosafe® 2 mL in 1/100 mL Serological Pipette	SINGLEPACK
OOPW-SP05	Oosafe® 5 mL in 1/10 mL Serological Pipette	SINGLEPACK
OOPW-SP10	Oosafe® 10 mL in 1/10 mL Serological Pipette	SINGLEPACK
Products available upon request. Please contact our Sales department: sales@sparmed.dk		
OOPW-CW03	Oosafe® Center Well Dish with two Compartments	MULTIPACK
OOPW-ST05	Oosafe® 60mm Dish, Treated, Label Area Grip	MULTIPACK
OOPW-TF05	Oosafe® 35mm Dish, High Wall, Treated	MULTIPACK

All plasticware products are packed in a secure white cardboard box.



facebook.com/SparMED



twitter.com/SparMED



linkedin.com/company/sparmed-aps

Oosafe® IVF TESTED

# PLASTIC WARE

FOR IVD USE



SparMED ApS  
Address: Ryttermarken 2,  
3520 Farum, DENMARK  
Phone number: +45 39-40-25-03



MK\_4.3.B2 rev. 01.06.20